Reliability and Mobility Meet Simplicity: MEDRAD® Salient

Clear Direction.  From Diagnosis to Care.

MEDRAD® Salient
Contrast Injection System

Clear Direction: MEDRAD® Salient Dual

The MEDRAD® Salient Dual Contrast Injection System has been engineered to combine innovation and mobility in order to bring simplicity and reliability to your suite.

Simplicity

- Intuitive, icon-driven interface
- Flexible workflow and software architecture
- Single-action syringe loading
- Easy installation
- Battery operation for total mobility

Reliability

- Allows dependable, reproducible flow rates and pressures that can consistently deliver injection protocols
- Provides a maximum pressure up to 300 psi and up to 10 mL/s
- Enables accurate and repeatable contrast administration with power injection¹⁻³

Benefits of Saline Flush

Reduce image artifacts

Using a saline flush helps minimize the streak artifacts that often occur as the superior vena cava empties into the right heart.⁴

Improve and increase peak contrast enhancement

Using a saline flush pushes the contrast through the vascular system.⁴

Controlled timing of contrast bolus

Using a saline flush creates a tight contrast bolus that enables precise timing of contrast in the vascular system, which is essential for multi-slice CT scanning.⁴

Optimize contrast usage

Using a saline flush can prevent pooling of contrast in veins and the contrast connection line.⁴

Versatility in Your Suite

The straightforward design of Salient Dual offers flexibility for your daily operations.

- Cordless configuration:
Removes tripping hazards
- Compact design:
Enhances mobility
Enables easy transport and storage
- Various positioning options:
Accommodates a range of spaces
- Customized protocols:
Implement customized protocols – up to six phases
Saves up to 20 customized injection protocols

Additional Tools for Enhanced Efficiency

Salient Dual can also support the following accessories and features:

- Wireless control room monitor
- Available scanner interface: ISI 700
- Protocol Assistant Tool (PAT)

Salient Workstation with the Protocol Assistant Tool (PAT)

Bayer Is Your Partner

MEDRAD® injectors can support you in achieving diagnostic image quality for your patients and workflow efficiency in your procedures.





Example injection protocol of fictional patient as shown on Control Room Monitor Screen

- The Protocol Assistant Tool (PAT) assists you to deliver customized contrast media protocols to your patients for many common CT examinations.
- Ultravist® (iopromide) works seamlessly with the Salient PAT software.
- Our Service and Applications Specialists provide quality service and solutions to meet your clinical, financial, and operational needs.
- To learn more about MEDRAD® Salient Dual and other Bayer products, please speak to your local representative.

Main Features and Benefits

Dimensions

MEDRAD® Salient D Technical Specifications

Functional

Maximum pressure	300 psi
Flow rate	0.1 to 10 mL/s
Maximum volume	190 mL
Multi-phase protocols	Up to six phases. Storage for 20 different named protocols
Manual fill rate	Up to 10 mL/s
Auto fill rate	Up to 10 mL/s

Controls

Base	Line power, IEC inlet, on/off switch
Head	5.7" LCD touch screen, arm button, on/off soft switch
Infrared remote control	Inject start/stop

Mechanical

Weight	25.5 kg (Dual)
Height	1315 mm
Floor area	500 mm x 500 mm
No. of wheels	4 casters, all individually locking

Electrical

Line voltage	100 – 230 VAC
Line frequency	50/60 Hz
Battery operation	Dual batteries, sealed lead acid
Mains power operation	Simultaneous operation and battery charging
Contrast heat maintainer (one (1) included)	37°C +/- 4°C, 5W

Ordering Details

PART NAME	PART NO.
Salient Single Injector	DC009S
Salient Single Injector With Remote Control Unit	DC009SW
Salient Dual Injector	DC009D
Salient Dual Injector With Remote Control Unit	DC009DW
Scanner Interface Box for ISI 700 (optional)	ISI2
Connector Tubing for Dual (box of 50)	ZY5152
Connector Tubing Single (box of 50)	ZY5151
Syringe With Quick Fill Tube (box of 50)	ZY6322
Syringe With Fill Spike (box of 50)	ZY6323
Syringe With Fill Spike, Single LPCT (box of 50)	ZY6324
Syringe With Quick Fill Tube, Single LPCT (box of 50)	ZY6325
Heat Maintainer (one (1) included)	DC022
Operation Manual – Multilanguage CD (one (1) included)	MN040015
*Single-Patient Line (250 cm)	SP0003
*Transfer Set and Connector Tubing	MP2002

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Clear Direction.  From Diagnosis to Care.

The patient data that appear in this document are fictitious protected health information (PHI) and/or actual PHI from which all personally identifiable information (PII) has been removed or otherwise anonymized. No personally identifiable information is shown.

Bayer, the Bayer Cross, MEDRAD Salient, MEDRAD, and Salient are trademarks owned by and/or registered to Bayer in the US and/or other countries. Other trademarks and company names mentioned herein are properties of their respective owners and are used herein solely for informational purposes. No relationship or endorsement should be inferred or implied.

© 2023 Bayer. This material may not be reproduced, displayed, modified, or distributed without the express prior written consent of Bayer.

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Bayer A/S
Imaxeon Pty Ltd
Bayer Medical Care B.V.